


How to Install the Citrix Web Plugin 12.3 for IE9


1. You will need local Administrator rights to perform this uninstall and reinstall of the Citrix client. Contact your local IT if you do not have local administrator access on your workstation.
2. Uninstall the current version of Citrix.
3. After Citrix has finished uninstalling the old Citrix client, go to www.7medical.com.
4. From the 7 Medical website select **“Support”**




5. Then go to the **“Download Center”**.




- From the Download Center select the **“Download”** button for **Citrix 12.3 Installer**.




- Select **“Run”** to start the install of the 12.3 Citrix client.



- Select **“Yes”** to allow access if the Window security window opens.
- Wait while the install is in progress.



- After the install is finished, select **“OK”**



- Launch Citrix.
- Login to the Citrix portal as normal, and launch your hosted application.