


X-Connect: One for All


The only integration engine you'll ever need.

Finally, you can manage all your interface needs using one multi-purpose integration engine, X-Connect™. Save time integrating among EMRs, billing systems, health information exchanges (HIEs), devices and more. X-Connect is the fastest, most efficient way to create, test, monitor and maintain interfaces, no matter the protocol or standard.

Health Data Repository Input Queue



Health Data Repository Output Queue



Any Data Format

HL7 v2, HL7 v3
FHIR
CDA
CCD
X12
NCPDP
DICOM
EDI
XML
Delimited Text

Any Protocol

LLP
HTTP/HTTPS
FTP/FTPS/SFTP
Web Services (SOAP, RESTful)
Database
File System
TCP/IP

WWW.7MEDICAL.COM

Key Benefits of X-Connect

There's no reason all your systems can't talk to each other

Too many disparate systems and too much duplication of data, makes for one big data headache. We feel your pain, and give you X-Connect. The all-in-one integration tool that solves your integration challenges, and maximizes data sharing, accuracy and workflow efficiency.

Simple Management Dashboard

Control all your interfaces and access stats from one centralized management dashboard. Manage workflows, set monitoring alerts, view real-time metrics.

Real-time Monitoring & Alerts

Quickly identify issues and exceptions so you can find and fix errors on the spot. Set alerts to escalate problems and configure workflows to assign troubleshooting tasks.

Easily Manage Users & Permissions

Quickly setup users and grant access permissions through a sophisticated user management module with advanced security hierarchy and rights options.

Efficient Mapping Tool

Easy to use tools for mapping data between systems regardless of the source, destination or data format. Reconfigure data mappings on the fly as requirements change.

Any System, Anywhere, Anytime

X-Connect exchanges data between any system (one-to-one, one-to-many and many-to-many). Access and manage interfaces from any location with Internet access. Full access to dashboards, logs and scripts from any Web-enabled device.

Around the Clock Service

A managed cloud service, X-Connect is backed by 24/7 support from healthcare workflow and technology experts to provide peace of mind.

Streamline your integrations. Get started today.

800.440.7119 | info@7medical.com | www.7medical.com

